Cochlear Vistafix®
Sterile products: Implants, abutments, and cover screws. Instructions for use

GB

Intended use: These products are used for craniofacial rehabilitation with bone anchored prostheses using the Vistafix System and should only be used by qualified professionals trained in the method ad modum Brånemark. Surgical literature and hands-on workshops are available from Cochlear. Go to www.cochlear.com for information.

Cochlear accepts no responsibility for any adverse outcomes if used with products not recommended by Cochlear. Close cooperation in an interdisciplinary team is essential for a successful outcome.

Indications: Craniofacial defects. Patients should have sufficient bone quality and quantity to support successful implant placement.

Contraindications:
a. Insufficient bone quality or quantity.
b. Patients (alone or with aid from others) that are unable to maintain and clean the skin around the abutment.

Children can be fitted, but treatment at an early age is dependent on the surgeon's evaluation regarding adequate bone depth. In the United States and Canada, the placement of a bone-anchored implant is contraindicated in children below the age of 5.

Possible complications: A 100 % rate of success cannot be guaranteed. Potential causes for failure of osseointegration include lack of adequate bone quantity/quality, trauma, infection, generalised diseases and surgical complications. A pre-operative evaluation determines possible risks from the implant procedure and healing capacities of the bone and/or associated soft tissue. Patients should be instructed to contact a physician in case of significant trauma, infection, skin irritation, pain, discomfort or implant loss. If the skin has been reduced and grows back to its original thickness, skin reduction surgery may be indicated.

Precautions: Before an MRI is performed on a patient with a Baha implant, read the MRI information supplied with this product. Only Cochlear products should be used with each other. However, combining products that are not dimensioned for mating can lead to mechanical failure, damage to tissue or to compromised results. Whenever the abutment needs to be adjusted, use appropriate wrench and counter torque wrench to minimise load on the implant.

Packaging and handling: All single use products are delivered sterile (sterilised using irradiation). Do not use products with damaged packaging or after their expiration date. Due to contamination and effectivity risks, do not re-sterilise or reuse these single use products.

Warranty: We, Cochlear, warrant to you, the consumer of the Product, that:
a. each Product is of merchantable quality;
b. each Product is reasonably fit for the purpose or purposes for which it is supplied by us; and
c. each Product will be free from defects in design, workmanship and materials for the Warranty Period. This warranty excludes liability for the defects or damage arising from, associated with, or related to the use of this Product with any non-Cochlear processing unit and/or any non-Cochlear implant. Warranty period for implants and abutments: 2 years starting on the date of surgery.

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