Cochlear  
Inform  
obrazo  
magne

## Cochlear™ Vistafix® MRI Safety Information

**The prosthesis, the bar construction, any fixation magnets, magnabutments or magnacaps attached to the implants must be removed before entering a room where an MRI scanner is located.**

Non-clinical testing has demonstrated that the Implant and Abutment [VXI300 & VXA300] are MR Conditional at 1.5 and 3.0 Tesla. It can be scanned safely under the following conditions.

- Static magnetic field of 1.5 Tesla and 3 Tesla only
- Maximum spatial gradient field of 3,000 Gauss/cm (30 T/m)
- Maximum MR System reported whole body averaged specific absorption rate (SAR) of 2 W/kg or maximum head averaged SAR of 3.2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the implant and abutment are expected to produce a maximum temperature rise of 1.1°C after 15 minutes of continuous scanning.

In non-clinical testing with the implant and abutment in place, the image artefact caused by the device extends approximately 1.8 cm from the implant and abutment when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. The artefact is reduced to 1.2 cm from the implant when the abutment is removed.

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